






































ROTORS AND ACCESSORIES







illustrated with carriers 5628
and lids 5629
with bio-containment⁵⁾
and removal frame 4626

Swing-out rotor Cat. No. (without carriers) 5624

capacity in ml	1.1–1.4	2.6–3.4	2.7–3	4–5.5	4.5–5	4.9	7.5–8.5	9–10	1.6–5	4–7	4–7	8.5–10
Ø x L in mm	8x66	13x65	11x66	15x75	11x92	13x90	15x92	16x92	13x75	16x75	13x100	16x100
Cat. No.	blood collection/urine tubes											
(adaptable only without lid 5629) adapting frame Cat. No. 5220-A  carrier Cat. No. 5628 ⁸⁾	 + 5220-A 	 + 5220-A 	 + 5220-A 	 + 5220-A 	 + 5220-A 	 + 5220-A 	 + 5220-A 	 + 5220-A 	 + 5220-A 	 + 5220-A 	 + 5220-A 	 + 5220-A 
Cat. No. ⁹⁾	5267	5268	5227	5264	5227	5268	5264	5258	5268	5264	5268	5258
boring Ø x L in mm	11 x 44	13.5 x 52	12.5 x 42	16 x 50	12.5 x 42	13.5 x 52	16 x 50	17.6 x 86	13.5 x 52	16 x 50	13.5 x 52	17.6 x 86
tubes per rotor	160	96	160	96	160	96	96	88	96		88	
max. RCF ²⁾	4,069	4,164	4,116		4,164		4,116	4,093	4,164	4,116	4,164	4,093
radius in mm	172	176	174		176		174	173	176	174	176	173
run-up in sec	90											
run-down in sec, braked	95											
temperature in °C ¹⁾	+4											

capacity in ml	15	50	50	250	two-tier frame for cyto accessories
Ø x L in mm	17x120	29x115	29x115	61x122	
Cat. No.	0509	-	0513	5127⁷⁾	cyto chambers¹⁰⁾
(adaptable only without lid 5629) adapting frame Cat. No. 5220-A  carrier Cat. No. 5628 ⁸⁾					
				+ 5220-A + 6319 	+ 5220-A 
Cat. No. ⁹⁾	6337-A	6339-A	6338-A	5263-A⁴⁾	5280
boring Ø x L in mm	17x90	30x87	30x90	66x87	-
tubes per rotor	56	24	24	4	16
max. RCF ²⁾ upper/lower	4,258	4,187	4,258	4,187	2,744 / 4,069
radius in mm oben/unten	180	177	180	177	116 / 172
run-up in sec	90				
run-down in sec, braked	95				
temperature in °C ¹⁾	+4				

D x W x H in mm	86 x 128 x 15	86 x 128 x 17.5	86 x 128 x 22
capacity in ml	-		
Cat. No.	MTP		CP
lid with bio-containment ⁵⁾ Cat. No. 5629  carrier Cat. No. 5628 ¹¹⁾			
Cat. No.	4626		
boring Ø x L in mm	-		
plates/strips pro rotor	24	20	16
max. RCF ²⁾	4,211		
radius in mm	max. 178		
run-up in sec	90		
run-down in sec, braked	95		
temperature in °C ¹⁾	+4		

⁷⁾ At temperatures of over + 40 °C and / or when not filled to capacity, bottles may warp during centrifugation.

⁸⁾ Please make sure not to exceed the max. permissible loading weight of 800g of carrier 5628. The loading weight is the sum of the weight of the respective adapters, of frame 5220-A (where required), the tubes, and their contents.

⁹⁾ With the exception of adapters 5263-A, 6337-A, 6338-A and 6339-A, two of these adapters need to be placed into carrier 5628.

¹⁰⁾ Please ask for our detailed cyto brochure.