Thermo Scientific Heraeus Megafuge Centrifuge Series





Specifications	Heraeus Megafuge 16/16R	Heraeus Megafuge 40/40R
Swinging Bucket Rotors		
Max Capacity	4 x 400 ml (with TX-400 rotor)	4 x 750 ml (with TX-750 and BIOLiner rotors)
Max Speed	5500 rpm (with TX-200 rotor)	5300 rpm (with BIOShield 1000A rotor)
Max RCF	5580 x g (with TX-200 rotor)	5590 x g (with BIOShield 1000A rotor)
Fixed Angle Rotors		
Max Capacity	6 x 100 ml (with F15-6x100 rotor)	6 x 100 ml (with F15-6x100 rotor)
Max Speed	15200 rpm (with 30x2 ml rotor)	15200 rpm (with 48x2 ml rotor)
Max RCF	25830 x g (with 30x2 ml rotor)	25314 x g (with 48x2 ml rotor)
Control System	Microprocessor	
Drive System	Direct, brushless induction low profile motor	
Rotor Locking System	Auto-Lock III	
Imbalance Detection System	SMARTSpin	
Programs	6, direct access	
Temperature Set Range*	-10 °C to +40 °C	
Pre-Cooling Function*	Yes, with direct button	
Refrigeration System*	CFC free	
Pulse (Short) Run	Yes	
Acceleration Rates	9	
Deceleration Rates	10	
Max Timer Range	9h, 99min + continuous	
Sound Level (dBA), with swinging bucket rotor at top speed	< 61/<55* dB (TX-400)	< 61/< 57* (TX-750)
Other Features	Certified ClickSeal biocontainment system, language selection, motorized lid latch	
Max Power Consumption (W) 120V	850/1200*	1400
230V	1010/1400*	1700/1950*
Max Heat Output (BTU/h) 120V	2900/4096*	4776
230V	3447/4776*	5800/6653*
Dimensions		
Height (lid open)	87 cm/ 34.3 in	90 cm/ 35.4 in
Height (lid closed)	36 cm/ 14.2 in	36 cm/ 14.2 in
Width	44 cm/ 17.3 in	62.3 cm/ 24.5 in
	62.3 cm*/ 24.5 in*	74.5 cm*/ 29.3 in*
Depth	60.5 cm/ 23.8 in	67 cm/ 26.4 in
Weight	57.5 kg/ 126.8 lb 91.5 kg*/ 201.7 lb*	86 kg/ 189.6 lb 116 kg*/255.7 lb*
Technical Standards		1 230 V only: EN 292, EN 61326, EN 55011B
Certifications	UL listed / CSA certified / CE marked / IVD compliant / Certified Biosafety	
Ordering Information	Heraeus Megafuge 16/16R	Heraeus Megafuge 40/40R
Heraeus Megafuge, ventilated, 230 V, 50/60 Hz	75004230	75004503
Heraeus Megafuge, ventilated, 120 V, 60 Hz	75004231	75004504
Heraeus Megafuge, refrigerated, 230 V, 50/60 Hz	75004251	75004518
Heraeus Megafuge, refrigerated, 120 V, 60 Hz	75004271	75004519

^{*}Refrigerated models









The quality of the procedures used to design, manufacture, sell and service the Heraeus Megafuge 16 and Megafuge 40 centrifuge series is certified and maintained in accordance to the **ISO 13485** international standard. Our effort to continuously improve our formalized business processes is certified and regularly audited in conformance with **ISO 9000** family of Quality Standards.

© 2009 Thermo Fisher Scientific Inc. All rights reserved. Corning and Corex are registered trademarks of Corning Inc. Kimble is a registered trademark of Kimble Glass. BD Falcon is a registered trademark of BD Biosciences. Sarstedt is a register trademark of Sarstedt, Inc. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America: USA/Canada +1 866 984 3766 (866-9-THERMO)

www.thermo.com/centrifuge

Ruth Allerta: 034/Callada +1 809 384 5706 (809-71EChNO)
Europe: Austria +43 1 801 40 0, Belgium +32 2 482 30 30, France +33 2 2803 2180, Germany national toll free 08001-536 376,
Germany international +49 6184 90 6940, Italy +39 02 02 95059 434-254-375, Netherlands +31 76 571 4440, Nordic/Baltic countries +358 9 329 100,
Russia/CIS +7 (812) 703 42 15, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220,

 $Other\ Asian\ countries\ +852\ 2885\ 4613\ \textbf{Countries}\ \textbf{not}\ \textbf{listed:}\ +49\ 6184\ 90\ 6940\ or\ +33\ 2\ 2803\ 2180$

